




# PRESSALIT®

<b>Product</b> <b>Item number</b>	CareSystem shower chair 410 with aperture, manually height adjustable R7167000
<b>Colour</b>          <b>Product description</b>	<div>  </div> White <p>Shower chair 410 with aperture, with back rest, "soft feel", seat height adjustable 340 mm.  Excl. mounting kit.  Discontinued as per 31/12-2015</p> <p>Shower chair 410 with aperture, with back rest, "soft feel", seat height adjustable 340 mm.  Excl. mounting kit.  Discontinued as per 31/12-2015</p>
<b>Adjustability</b>  <b>Product capacity</b>	Height adjustable  Max. vertical load: 150 kg.
<b>Materials</b>          <b>Ergonomy</b>	Basic material: 20 micrometre anodized aluminium  Seat/backrest: Polyurethane foam with marine ply core  Seat/backrest surface: Gum lac  Lever: Stainless steel w/nylon knob  Front plate/end covers, track: Polystyrene plastic (PS)  Soft and slip-resistant material with a warm feel.  Ergonomically designed seat with aperture that makes intimate hygiene easier.



# PRESSALIT®

<b>Product</b> <b>Item number</b>	CareSystem shower chair 410 with aperture, manually height adjustable R7167000
<b>Features</b>	<p>The seat can fold up when not in use or to make more space.</p> <p>Height adjustable 340 mm - the height of the seat and the backrest can be adjusted separately.</p>
<b>Cleaning</b>	Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®