




PRESSALIT®

Product Item number	CareSystem shower chair 410, manually height and sideways adjustable R7566000
Colour Product description	<div>  </div> White Shower chair 410 with back rest and arm supports, "soft feel", seat height adjustable 330 mm, for mounting on horizontal track. Discontinued as per 31/12-2015 Shower chair 410 with back rest and arm supports, "soft feel", seat height adjustable 330 mm, for mounting on horizontal track. Discontinued as per 31/12-2015
Adjustability Product capacity	Height and sideways adjustable Max. vertical load: 150 kg.
Materials Ergonomy	Basic material: 20 micrometre anodized aluminium Seat/backrest: Polyurethane foam with marine ply core Seat/backrest surface: Gum lac Armrest surface: Two-component plastic (PVC) End cap, armrest: Arnite Lever: Stainless steel w/nylon knob Front plate/end covers, track: Polystyrene plastic (PS) The armrests' oval shape gives an easy and comfortable surface to lean and push against. Soft and slip-resistant material that has a warm feel. Ergonomically designed seat.



PRESSALIT®

Product Item number	CareSystem shower chair 410, manually height and sideways adjustable R7566000
Features	<p>The seat and armrests can fold up when not in use or to make more space.</p> <p>Height adjustable 330 mm - the height of the seat and the backrest can be adjusted separately.</p> <p>Can be pushed along the track.</p>
Cleaning	Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®