




# PRESSALIT®

<b>Product</b> <b>Item number</b>	CareSystem front and back support, 700 - 770 mm RT802000
<b>Colour</b>  <b>Product description</b>	White  <p>Front and back support for shower chair, adjustable A: 700 mm - 770 mm. For support arms with centre distance of 700 mm.</p> <p>Front and back support for shower chair, adjustable A: 700 mm - 770 mm. For support arms with centre distance of 700 mm.</p> 
<b>Product capacity</b>	Max. vertical load: 70 kg.
<b>Materials</b>	Two-component plastic (PVC) Strap: stainless steel
<b>Ergonomy</b>	Stabilises and corrects a user with reduced body balance.
<b>Features</b>	Adjustable 700-770 mm.
<b>Cleaning</b>	Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.



# PRESSALIT®

<b>Product</b>	CareSystem front and back support, 700 - 770 mm
<b>Item number</b>	RT802000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®