




PRESSALIT®

Product Item number	CareSystem backrest, sideways adjustable R6513000
Colour Product description	<div> <div>White</div>  </div> <p>Discontinued 31 Dec 2015. Back rest for shower chair for toilet without cistern. For mounting on horizontal track. (recommended for toilet pans with a total distance of 700 mm from wall to front)</p> <p>Back rest for shower chair for toilet without cistern. For mounting on horizontal track. (recommended for toilet pans with a total distance of 700 mm from wall to front)</p>
Adjustability Product capacity	Sideways adjustable Max. vertical load: 100 kg.
Materials Ergonomy Features	<p>Basic material: 20 micrometre anodized aluminium</p> <p>Backrest: Polyurethane foam with marine ply core</p> <p>Backrest surface: Gum lac</p> <p>Lever: Stainless steel w/nylon knob</p> <p>Soft and slip-resistant material with a warm feel.</p> <p>For toilet without cistern.</p> <p>Recommended when the toilet's projection = 700 mm.</p> <p>Can be pushed along the track.</p>
Cleaning	<p>Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.</p>



PRESSALIT®

Product	CareSystem backrest, sideways adjustable
Item number	R6513000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®