




PRESSALIT®

Product Item number	CareSystem front and back support, 540 - 620 mm RT801000
Colour Product description	White Front and back support for shower chair, adjustable A: 550 mm - 620 mm. For shower chair 410 and support arms. Front and back support for shower chair, adjustable A: 550 mm - 620 mm. For shower chair 410 and support arms.
Product capacity	Max. vertical load: 70 kg.
Materials	Two-component plastic (PVC) Strap: stainless steel
Ergonomy	Stabilises and corrects a user with reduced body balance.
Features	Adjustable 550-620 mm.
Cleaning	Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.





PRESSALIT®

Product	CareSystem front and back support, 540 - 620 mm
Item number	RT801000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®